

OsmoTECH® PRO Multi-Sample Micro-Osmometer

Designed to meet the most demanding needs in your Bioprocess, the OsmoTECH PRO delivers the highest level of data management, accuracy and precision, and ease-of-use available on the market.

- Drive **unparalleled operational efficiency** with an intuitive interface and 20-position walkaway testing design.
- Easily integrate data into your established workflow with **unsurpassed data management capabilities**.
- Evolve with the changing needs of your organization with **highly configurable settings**.
- Stay in compliance with the growing number of regulations with the system's **robust data integrity features**.



Dependable. Trusted. Designed for Biotech.



Flexible data management options to fit your specific workflow

- ✓ **OPC-UA:** Easily export your results into a data management system and batch records
- ✓ **PDF or .CSV export:** Quickly save data to network database or USB
- ✓ **Database backup:** Protect your data with automatic or manual backups
- ✓ **Synchronize with network time:** Eliminate the need to manually manage date and time
- ✓ **Embedded Web server:** Remotely view and manage data for additional walkaway time




Security features support 21 CFR Part 11 and EU Annex 11 compliance

- ✓ **3 levels of user access:** Set up to reflect the hierarchy of your team
- ✓ **Active Directory / LDAP:** Align with your IT network to minimize manual input
- ✓ **Configurable password strength:** Ability to match internal security requirements
- ✓ **Unlimited results and events storage:** Preserve your audit trail
- ✓ **Pre-expired passwords:** Force password change on first login to ensure maximum security
- ✓ **Require reason for settings changes:** Maintain a thorough audit trail
- ✓ **Comments field:** Add unique comments to test results for complete audit trail
- ✓ **Automatic logout timer:** Keep your data safe during idling
- ✓ **Results database on/off:** Decide whether results are stored on the system

Learn more: acompanies.com/OsmoTECHPRO

By implementing OsmoTECH PRO into your workflow, you can:



Parts and supplies

Part number	Product description
222825	Disposable sample tubes
222840	Probe wiper rings
800097	Pipette tips
202850	Swab cleaner kit
3MA005	50 mOsm/kg Calibration Standard
3MA085	850 mOsm/kg Calibration Standard
3MA200	2000 mOsm/kg Calibration Standard (optional)
3MA029	Clinitrol™ 290 Reference Solution
135022	Dot matrix printer
AN2TP5	Printer paper

Specifications

Sample type	Aqueous solution
Sample volume	30 ± 1 µL
Test time	90 seconds
Sample capacity	20 samples
Range	0 to 2000 mOsm/kg H ₂ O

Not for patient diagnostic use



ADVANCED
INSTRUMENTS

Two Technology Way / 781-320-9000
Norwood, Massachusetts 02062, USA
800-225-4034 Fax: 781-320-8181
aicompanies.com

©2019 Advanced Instruments. Clinitrol and OsmoTECH are trademarks of Advanced Instruments. All other trademarks are the property of their respective companies. Advanced Instruments certifies that the technical features needed for 21 CFR Part 11 and EU Annex 11 compliance are built into OsmoTECH PRO. It is your responsibility to implement the necessary controls in your laboratory to comply with 21 CFR Part 11 and EU Annex 11 requirements.

PCN00941

MP00065 Rev1

Request a quote today: aicompanies.com/OsmoTECHPRO-quote